


## Cells on chips

Mol3014, Øyvind Halaas, IKM

[www.ntnu.no](http://www.ntnu.no)


### Course outline : nanotechnology for biomedicine : analysis



### Learning goals

- Describe shortcoming of current methods of cell biology
- Elaborate on why single cell analysis is necessary
- Provide examples on cell standardization by micropatterning
- Provide examples of cell traps
- Show possibilities for cell-pairing
- Provide examples on cell sorting in microdevices
- Describe devices for migration in channels and on patterns

# We are cells



## Cells

The minimal unit of life

About 200 standard cell types

Neurons and immune cells are highly dynamic

Immune cells are highly diverse (many millions genotypes)

Cancer cells are changing and non-clonal (metastatic cancer cell very different from initial cancer stem cell)

## Cell biology

The study of cells provide information about

Gene function (research, about 50% genes unknown function)

Development (embryology, tissue repair)

Stem cells (growth requirements, growth patterns)

Immune cells (vaccines, autoimmunity, cancer, allergies)

Neurons (cell signalling and networks)

Cancer cells (growth, aggression, metastates, drug targets)

Drug response (activity testing)

## What can we learn from cell studies

### **Birth (proliferation)**

cancer, tissue regeneration, immunology

### **Food (metabolism, pinocytosis, endocytosis, phagocytosis)**

disease, cancer, immunology

### **Action (motility, movement)**

immunology, cancer metastasis, tissue (re)generation

### **Looks (morphology)**

gene function, development, neuronal network

### **Social networking (signalling)**

gene/protein level alteration, internal body communication

### **Death (apoptosis)**

cancer, immunology, development, disease

## **Readouts**

### **Quantitation and Identification of**

Metabolites (intracellular, extracellular, enzyme inhibitors)

Nucleotide (DNA, mRNA, miRNA)

Protein (eg cytokines, induced/repressed proteins)

Cell (cell subtype, function, morphology, quality)

Interactions between them

Location, timing, levels, networks


## **Methods**

Imaging (morphology, localization)

Biochemical (identity and quantity)

Informatical (system)

## What to analyze and how



Measurement of single-cell dynamics  
David G. Spiller, Christopher D. Wood, David A. Rand & Michael R. H. White  
Nature 465, 736–745 (10 June 2010) doi:10.1038/nature09232

## Motivation for new approaches:

**Today's bulk approaches have severe practical limitations and do not mimic tissues satisfactorily**

**Cells are unsynchronous, heterogeneous and depend on ECM and cell:cell contact**


**Today we have no spatial control and poor temporal control**

**Single cell studies is today uncoupled from single cell analysis (standard microscopy exception)**


**Solutions:**

Single cell analysis  
Standardization of cells  
Microfluidics  
Microdroplets

Cells and even cell lines are heterogenous, this is averaged by bulk research, but represent problems for rare events and unsynchronized cells.




DOI: 10.1016/j.tibtech.2010.03.002



## How to affect cell fate


Microtechnologies offers the opportunity for highly localized cell perturbations




Biotechnol J. 2013 Feb;8(2):167-79.  
doi: 10.1002/biot.201200149. Epub  
2012 Aug 14.  
Arrayed cellular environments for  
stem cells and regenerative medicine.  
Titmarsh DM, Chen H, Wolvetang EJ,  
Cooper-White JJ.

Single cell studies of today

A confocal only allow focus on a single layer determined by pinhole






Cells are annoyingly hetergeneous



Ø Halaas. 2009




**Cell imaging demands experienced researchers.  
Manuscript referees demand quantitative proof**



Question: How often does the red drug delivery system meet the green intracellular bacteria?

Optical imaging essential, but new tools are needed to become objective

## Cells are dynamic structures



seconds

minutes

hours

Changes on many time scales



## Microscopes

Good at high content (morphological detail) analysis  
Good at high resolution (spatial) analysis  
Good at time-resolved analysis

Expensive  
Bad at analysing many cells  
Information only for 3-4 different molecules


Single cell tools of today, the flow cytometer

**A flow cytometric cell sorter  
(1990)**




What you get when crossing an ink-jet printer with Coulter counter (particle size analyzer), Palo Alto 1965

## A flow cytometric single cell sorter



## Flow cytometry can detect rare cell types

One can analyze up to 8 protein levels in each cell subset



## Flow cytometer

- Good at analysing many cells (millions)
- Fast instrument (20.000 cells/s)
- Can analyze up to 8 fluorescent parameters + size and optical density
- Bad at dynamic analysis (is mostly end-point analysis)
- Expensive
- No morphological information

## Next generation cell biology: Cells on chips




Figure 2





Illustration of "small-volume effects" in microfluidic cell culture devices. Smaller culture media volume for a given cell results in faster consumption of nutrients and increased concentration of metabolites or secreted molecules, similar to tightly packed tissues.

Lab Chip. 2013 Aug 27;13(19):3789-802. doi: 10.1039/c3lc50550b.  
**Cell chips as new tools for cell biology - results, perspectives and opportunities.**  
 Primiceri E, Chiriaco MS, Rinaldi R, Maruccio G.

What is missing from conventional cell cultures?



Spatial and temporal control of cell environment  
at the scale of the cell (0.1 to 100  $\mu\text{m}$ , ms to days)

**a) Control of growth medium** (microfluidic):

- chemostat, thermostat
- drug delivery (serial dilutions, oscillations, etc...)
- gradients

**b) Surface control** (micropatterning, micro-fabrication):

- adhesion (stamping of adhesion molecules, contact guidance)
- cell shape
- cell positioning (for cell/cell communication)

**c) Substrate mechanical properties**

**d) Volume control for 3D confinement** (microchambers for cells or in-vitro assays)

**And all combinations...**

Table 1 Selected Studies on Stem Cells Using Lab-on-a-chip Devices

Stem Cell Application	Selected Studies
Control of Soluble Factors	Dynamic exposure to soluble drug Design for controlled-diffusive mixing Microfluidic based multi-injector method 'Microbioreactor' system for human embryonic stem cell differentiation Paracrine and autocrine signaling control
Control of ECM Interactions	Combining soluble and insoluble factors in microarray 3D microenvironment to study stem cell fate
Control of Cell-Cell Interactions	Design to control cell-cell interactions Design for cell pairing Design for co-culture of mES cells with other cell types Design for temporal control of cell-cell signaling Design for heterotypic cells with different shear sensitivities
Control of Mechanical Signals	Controlled elasticity geometry through PDMS patterns Control of elasticity using magnetic nanoposts Controlling nanoscale symmetry and disorder in cues Designs for surface pattern geometries Design for study of role of topography Control of 3D topography by electrospun nanowires Control of shear stress as a mechanical cue Shear control for long term culture
High Throughput Screening	High throughput integrated microfluidic systems 'High-content screening' with combinatorial stimuli Design to study high throughput temporal responses Design to screen effects of drugs in long term cultures
Novel Studies for Potential Stem Cell Research	Slow perfusion by osmotic pump for long term culture Design for reversible bonding of microfluidic device to cell culture containing cover slip by vacuum Design for reversible bonding of microfluidic device by magnetic force



Table 1 Important differences between conventional (monolayer) and microfluidic *in vitro* cell culture systems

Conventional	Microfluidics
No confinement (open wells).	<i>Cell microenvironment</i> Confinement (closed systems).
Limited level of <b>spatial control</b> (e.g. only single-well or trans-well systems).	High level of <b>spatial control</b> (e.g. compartmentalization for co-culture, 3-dimensionality and sub-cellular resolution).
No <b>fluid control</b> (only static or chaotic).	High level of control over fluids (e.g. laminar flow, perfusion, and temporal control over fluid exchange).
Limited possibilities for creating <b>physical stimuli</b> .	Various <b>physical stimuli</b> possible (e.g. stiffness, shear, compression).
Low temporal and spatial control over <b>chemical stimuli</b> (only bulk addition).	Possibility to create highly defined spatial and temporal <b>chemical stimuli</b> (e.g. soluble or surface gradients).
Established and characterized culture <b>substrate materials</b> .	Limited characterization of applied <b>substrate materials</b> and limited use of biologically characterized materials.
Compatible with conventional standardized <b>biological assays</b> .	<i>Biological read-out</i> Compatibility issues with conventional standardized <b>biological assays</b> .
Compatible with established <b>read-out equipment</b> .	Compatibility issues with established <b>read-out equipment</b> .
Comparable to large amount of data from <b>historical experiments</b> .	Low number of available <b>historical experiments</b> limits comparison.
Limited possibilities for <i>in situ</i> <b>read-out of biological processes</b> .	Possibility to integrate sensors and assays for <i>in situ</i> <b>read-out of biological processes</b> .
High <b>reagent</b> and biological (cell) <b>material use</b> in HTS setting.	<i>High-throughput screening (HTS)</i> Reduced <b>reagent</b> and biological (cell) <b>material use</b> in HTS setting.
Limited possibilities to <b>parallelize and integrate</b> assays.	Highly applicable to <b>parallelization and integration</b> of assays.
Compatible to conventional high-throughput (robotics) <b>equipment</b> .	Not compatible with conventional high-throughput (robotics) <b>equipment</b> .

# Single cell devices

- Watching single event dynamics
- Relate to earlier lectures
- 1) Standardize single cell morphology
- 2) Isolate cells
- 3) Trap cell
- 4) Cell pairing
- 5) Select cells
- 6) Cell migration
- 7) Multiplexing
- 8) Cocultures


Application	Challenges with traditional methods	Advantages of microfluidics
Single-cell RT-qPCR	Limited abundance of starting template Cost and throughput required for analysis of large numbers of cells and target genes	Concentration enhancement in small volumes Parallelization, automation, and economy of scale
Single-cell genomics	Amplification bias and sensitivity	Improved reaction bias and sensitivity in nL volumes and reduced contaminant DNA
Single-cell measurements of intracellular proteins	Isolating individual cells Movement of living cells	Integrated microfluidic cell sorting and processing Confinement of live cells in microfluidic structures
Single-cell measurements of secreted proteins	Low amount of signal Small amounts of secreted products from single cells Difficult to co-localize multiple cells in defined chemical environments	Integrated single-cell handling allows lysate analysis Concentration enhancement in small volumes Ease of confinement in droplets or microchambers
Signaling studies	Mostly limited to static conditions Inability to rapidly exchange conditions on suspension cells	Easy temporal stimulation Laminar flow and proper design enables cell sequestration
Live cell imaging	Difficulties of tracking cells through multiple frames	Confinement of clones facilitates cell tracking



## 1) Standardization of single cell organization




<http://cytoo.com/>

## Standardizing cell cultures




Nelson, C.M., Science, 2006

## 2) Single cell capture in microwells





2) Microwells with detection modes



Some cells are not adherent and dont want to be:  
Alternative to adherent cells: Hanging drop cultures in microfluidics





Figure 1: A cross-section of the operation of the microfluidic HD platform. (a) The microfluidic channel is filled and droplets are formed and hanging at the openings. (b) Cells are driven by the flow and docked in the HD. (c) Untrapped cells are rinsed off. (d) Cells are incubated and their extracellular solution can be modulated through the microchannel while gas can also be exchanged directly through the HD that is open to air.

Most cells are not in their natural state on a hard surface  
- Suspended single cell cultures




Nanoscribe, direct 3D writing

### 3) Cell traps

Alternative to  
printed biopatterns

Capture of non-  
adherent cells at  
specified locations  
exploiting fluid  
dynamics




### 3) Cell traps


Nat Methods. 2009 Feb;6(2):147-52. doi: 10.1038/nmeth.1290. Epub 2009 Jan 4.

**Microfluidic control of cell pairing and fusion.**

[Skelley AM](#), [Kirak O](#), [Suh H](#), [Jaenisch R](#), [Voldman J](#).




### 3) A resettable cell trap : how to get cells out of the chip after trapping





Ryan D. Sochol<sup>1\*</sup>, Kosuke Iwai<sup>1</sup>, Adrienne T. Higa<sup>1</sup>, Joanne C. Lo<sup>2</sup>, Eric Zhou<sup>1</sup>, Liana Lo<sup>1</sup>,  
Chris Luong<sup>1</sup>, Megan Dueck<sup>1</sup>, Song Li<sup>1</sup>, Luke P. Lee<sup>1</sup> and Liwei Lin<sup>1</sup>  
Berkeley Sensor and Actuator Center  
<sup>1</sup>University of California, Berkeley USA  
<sup>2</sup>Sandia National Laboratories, Livermore USA


### 3) Using microfluidics for trapped single cell injections




Microfluidic based single cell microinjection.  
Lab Chip. 2008 Aug;8(8):1258-61.  
Adamo A, Jensen KF.



### 3,4) Cell trapping and pairing example from hybridomas




Hybridomas  
Immune cells  
Stem cells  
Cancer cells  
Neurons



[http://www.nature.com/nmeth/journal/v6/n2/supinfo/nmeth.1290\\_S1.html](http://www.nature.com/nmeth/journal/v6/n2/supinfo/nmeth.1290_S1.html)

#### 4) Microfluidic cell pairing




A microfluidic array with cellular valving for single cell co-culture  
 Jean-Philippe Frimat, Marco Becker, Ya-Yu Chiang, Ulrich Marggraf, Dirk Janasek, Jan G. Hengstler, Joachim Franzke and Jonathan West  
 Lab Chip, 2011, 11

DOI: 10.1039/C0LC00172D


#### 5) Cell isolation

- Adhesion
- Size filtering
- Dielectrophoresis
- Fluorescence

# Circulating tumor cells



## Liquid biopsies, circulating tumor cells



CTCs are few and heterogenous, needs to be studied at single cell level

- 1) Enrichment
- 2) Single cell analysis

Lab Chip. 2013 Aug 28. [Epub ahead of print]  
**Advances and critical concerns with the microfluidic enrichments of circulating tumor cells.**

[Hyun KA, Jung HI.](#)






## 5) Serial size-filtering of CTC


[Detection and isolation of circulating tumor cells: Principles and methods.](#)

Esmailsabzali H, Beischlag TV, Cox ME, Parameswaran AM, Park EJ. Biotechnol Adv. 2013 Aug 30. doi:pii: S0734-9750(13)00150-X. 10.1016/j.biotechadv.2013.08.016.


<http://www.youtube.com/watch?v=xmgP8fBfcJU>



## 5) On chip filtering capture and sample prep for larger cells




## 5) Micro/nanotechnological options for CTC enrichment



## 6) Migration

- Stem cell for tissue engineering (in vitro)
- Stem cells for tissue regeneration (in vivo)
- Tissue development and remodelling
- Immunity
- Cancer metastasis and tissue invasion
- Other diseases

# 6) Migration/Chemotaxis




Lab Chip. 2013 Jul 7;13(13):2484-99. doi: 10.1039/c3lc50415h. Epub 2013 May 28.  
Recent developments in microfluidics-based chemotaxis studies.  
Wu J, Wu X, Lin F.


# 6) Soluble Gradients

Migration  
towards  
gradients

Constant flow  
or  
immobilize with  
photochemistry



## 6) Migration




Lab Chip. 2013 Aug 27;13(19):3789-802. doi:  
10.1039/c3lc50550b.

**Cell chips as new tools for cell biology -  
results, perspectives and opportunities.**



Primiceri E, Chiriaco MS, Rinaldi R, Maruccio G.

### 6) Migration in microchannels:

microchannels with a 4µm by 5µm section




Science. 2008 Dec 12;322(5908):1705-10.  
Regulation of dendritic cell migration by CD74, the MHC class II-associated invariant chain.  
Faure-André G, Vargas P, Yuseff M, Heuzé M, Diaz J, Lankar D, Steri V, Manry J, Hugues S, Vascotto F, Boulanger J, Raposo G, Bono MR, Roseblatt M, Piel M, Lennon-Duménil AM.



<http://www.sciencemag.org/cgi/content/full/sci;322/5908/1705/DC1>

Faure-André et al., Science, 2008



**Own research: Recreating infections on a chip**  
**Understanding recruitment toward tuberculosis granuloma**



Phagocytic cells stop moving when dealing with substrates

Mycobacteria are resilient.

Eating dead cells=cross presentation




Cell retention filters segregates compartments



Moving in and out of the artificial tissue

Forcing cells to make choices where we control input



6) Migration on tracks (not tunnels)

J Cell Biol. 2009 February  
23; 184(4): 481–490.  
doi:  
10.1083/jcb.200810041






Figure 2.  
1D topography functionally mimics 3D fibrillar matrix. (A) Rapid migration of NIH-3T3 fibroblasts through a 3D cell-derived matrix. (B) Cytoskeletal alignment (green, actin) along oriented FN fibers (red). (C) Fibroblasts plated on single ~1.5-μm lines mimic this 3D phenotype and migrate continuously. (D) Cell migration rates on 2D surfaces and 1D fibrillar lines at different FN ligand densities versus 3D cell-derived matrix. (E) HK migration over 2D, 3D matrix, and 1D fibrillar lines. (F) HK migration rates on substrates in E. Bars: (A, B, and E) 20 μm; and (C) 10 μm. \*, P < 0.01 versus 2D; †, P < 0.05 versus 0.5 and 1,000 FN. Error bars indicate SEM.

6) Also adhesive patterns in the right shape can induce motility






## Conclusions from extensive review 1D-2D-3D




## Vertical compartmentalization

A tool for looking at migration  
in blood-model






Chem Soc Rev. 2010 Mar;39(3):1036-48. doi: 10.1039/b909900j. Epub 2010 Feb 1.  
Fundamentals of microfluidic cell culture in controlled microenvironments.  
Young EW, Beebe DJ.

## Combining confinement and micro-fluidics: dual chambers for neuronal growth



Taylor et al., 2005


## Neurons on chips



## 7) Multiplexing cell cultures

Giving different mixtures of different biologicals to different cells,


Robotless- High throughput screening on a chip




**Scheme 1** Chip structure. A: Illustration for the three layered components of the chip: the top two are PDMS layers, and the bottom one is the cell culture substrate layer. B: Illustration for the structure of integrated chip. C: Illustration for the single reservoir unit of the chip.

The third dimension gaps cell culture with tissues (more next lecture «organs on chips»)

<http://www.ncbi.nlm.nih.gov/pubmed/22797912>




## 8) Options for in vitro co culture cell growth



## 8) New approaches for cell cultures

Curr Opin Biotechnol. 2011  
Oct;22(5):690-7. doi:  
10.1016/j.copbio.2011.07.00  
1. Epub 2011 Aug 5.  
Biomaterials meet  
microfluidics: building the  
next generation of artificial  
niches.  
Kobel S, Lutolf MP.



## General conclusions

A very large variety of techniques: how to choose the one you need?

**1. Cell type:** imposes the type of protein, and the type of protein binding

**2. Length of the experiment:** imposes the type of anti-adhesive molecule and the binding method

### 3. Main parameters to modulate

- Size of the patterns (bigger are easier, threshold is around few tens of microns)
- Cell culture substrate (glass is the best for optics and gold is the best for chemistry)
- Anti-adhesive backfilling (Silanes, PLL-PEG, Pluronic, BSA, etc...)
- Protein binding to the substrate (covalent, electrostatic or adsorption)

Some methods, like photolithography and lift-off work in almost all cases, but are very heavy for a biology lab, while others are very easy but can apply only to very specific cases

pensum

Analyst. 2016 Jan 21;141(2):504-24. doi: 10.1039/c5an00382b.

Microfluidics for research and applications in oncology.

Chaudhuri PK, Ebrahimi Warkiani M, Jing T, Kenry, Lim CT.

Overview of single-cell analyses: microdevices and applications  
Sara Lindström<sup>a</sup> and  
Helene Andersson-Svahn  
Lab Chip, 2010,10, 3363-3372

Microfluidics for single cell analysis.  
Yin H, Marshall D.  
Curr Opin Biotechnol. 2012 Feb;23(1):110-9. doi:  
10.1016/j.copbio.2011.11.002. Epub 2011 Nov 29. Review.

<http://www.youtube.com/watch?v=66Oc8fLKKIE>